

2- door ventilated under-bench safety cabinet - 2 compartments, including PVC equipments for acids

COMPLIANCE

ADVANTAGES

- Made in France
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Normalized labelling provided with the equipment
- Holes for the connection to the ventilation
- Allows the separation of different types of chemicals
- PVC protections for acid compartment (included)



TECHNICAL FEATURES

Model	AB56	Ventilation	Inlet located at the bottom of the cabinet, outlet on the back wall with flame arresters
Cabinet materials	10/10th epoxy painted steel structure	Openings diameter (mm)	100
Doors color	white RAL 9010	Number of compartments	2
Number of doors	2	Shelves in PVC	1
Doors type	swing doors	Containment sump	1
Hinges	right & left	Containment sump in PVC	1
Door closure	manual closing		50.00
Shipping dimensions HxWxD (mm)	870 x 1200 x 800	HS Code	94032080

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Connecting kit (Réf : KRC)
- Vertical recycling air box (delivered without filter) (Réf : H40+)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
AB56	670 x 1100 x 545	2x (645 x 540 x 475)	60		2	50	50.00