

2- door under-bench safety cabinet - 2 compartments, including PVC equipments for acids

COMPLIANCE

ADVANTAGES

- Made in France
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- PVC protections for acid compartment (included)
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Allows the separation of different types of chemicals



TECHNICAL FEATURES

Model	AB56SV	Ventilation	Inlet located at the bottom of the cabinet, outlet on the back wall with flame arresters
Cabinet materials	10/10th epoxy painted steel structure	Openings diameter (mm)	100
Doors color	white RAL 9010	Number of compartments	2
Number of doors	2	Shelves in PVC	1
Doors type	swing doors	Containment sump	1
Hinges	right & left	Containment sump in PVC	1
Door closure	manual closing		50.00
Shipping dimensions HxWxD (mm)	870 x 1200 x 800	HS Code	94032080

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Connecting kit (Réf : KRC)
- Vertical recycling air box (delivered without filter) (Réf : H40+)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
AB56SV	670 x 1100 x 545	2x (645 x 540 x 475)	60		2	50	50.00