



2-door large sheet steel safety cabinet

COMPLIANCE

ADVANTAGES

- Compliant to the Public Health Code
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Locking system for the other drawers in the retracted position when one drawer is already pull-out
- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- Containment sumps on sliding rails
- INRS recommendations compliant
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Possibility to fix on the wall



TECHNICAL FEATURES

Model	ALT305	Shipping dimensions HxWxD (mm)	2150 x 1200 x 800
Cabinet materials	8/10th epoxy painted steel structure	Ventilation	Both inlet and outlet
Doors color	grey RAL 7035	Openings diameter (mm)	100
Number of doors	2	Number of compartments	1
Doors type	swing doors	HS Code	102.00
Hinges	right & left		94032080
Door closure	manual closing		

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Containment sump for 2-door cabinet (Réf : ESALT305)
- Containment sump in PVC for 2-door cabinet (Réf : SEALT305)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
ALT305	1950 x 916 x 500	1718 x 795 x 430	300		0	50	102.00