

AMS100



Safety cabinet for flammable products - With CO2 extinguisher

COMPLIANCE

CNPP

ADVANTAGES

- Retention equipment
- CNPP compliant
- Storage of flammables

BUT ALSO...

- Sloping roof to avoid unauthorized storage
- Sloping roof to avoid unauthorized storage
- CNPP compliant
- Locking system
- Self-closing doors
- Big containment sump
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Equipped with a CO2 extinguisher
- Automatic release by thermic bulb (between 65 and 72°C)



TECHNICAL FEATURES

Model	AMS100	Shipping dimensions HxWxD (mm)	1700 x 1200 x 800
Compliance	CNPP	Ventilation	Both inlet and outlet
Cabinet materials	15/10th epoxy painted steel structure	Openings diameter (mm)	100
Doors color	yellow RAL 1023	Number of compartments	1
Number of doors	2	Containment sump	1
Doors type	swing doors	HS Code	70.00
Hinges	right & left		94032080
Door closure	self-closing		

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Additional shelf for AMS100 - AMS100SE (Réf : ESAMS1)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
AMS100	1500 x 800 x 400	1250 x 780 x 370	100		3	50	70.00