

## Safety cabinet for flammable products - Without extinguisher

### COMPLIANCE

CNPP

### ADVANTAGES

- Retention equipment
- CNPP compliant
- Storage of flammables

### BUT ALSO...

- Sloping roof to avoid unauthorized storage
- Sloping roof to avoid unauthorized storage
- Locking system
- CNPP compliant
- Self-closing doors
- Normalized labelling provided with the equipment
- Holes for the connection to the ventilation
- Big containment sump



### TECHNICAL FEATURES

Model	<b>AMS100SE</b>	Shipping dimensions HxWxD (mm)	<b>1700 x 1200 x 800</b>
Compliance	<b>CNPP</b>	Ventilation	<b>Both inlet and outlet</b>
Cabinet materials	<b>15/10th epoxy painted steel structure</b>	Openings diameter (mm)	<b>100</b>
Doors color	<b>yellow RAL 1023</b>	Number of compartments	<b>1</b>
Number of doors	<b>2</b>	Containment sump	<b>1</b>
Doors type	<b>swing doors</b>	HS Code	<b>70.00</b>
Hinges	<b>right &amp; left</b>		<b>94032080</b>
Door closure	<b>self-closing</b>		

### ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

### OPTIONS

- Additional shelf for AMS100 - AMS100SE (Réf : ESAMS1)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
AMS100SE	1500 x 800 x 400	1250 x 780 x 370	100		3	0	70.00