

2- glass door tall safety cabinet

COMPLIANCE

ADVANTAGES

- Made in France
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Locking system



TECHNICAL FEATURES

Model	ASDV	Shipping dimensions HxWxD (mm)	2025 x 1200 x 800
Cabinet materials	10/10th epoxy painted steel structure	Ventilation	Inlet located at the bottom of the cabinet, outlet located at the top
Doors color	white RAL 9010	Openings diameter (mm)	100
Number of doors	2	Number of compartments	1
Doors type	swing doors	Containment sump	1
Hinges	right & left		95.00
Door closure	manual closing	HS Code	94032080

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Retention shelf for ASDP - ASDD - ASDV (Réf : ESASDD)
- Containment sump in PVC (Réf : SBPVC)
- Shelf protection in PVC (Réf : SEPVC)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
ASDV	1825 x 1100 x 545	1800 x 1090 x 475	300		3	0	95.00