

2- glass door tall safety cabinet - 2 compartments

COMPLIANCE

ADVANTAGES

- Made in France
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Glass doors for an immediate view of the products
- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Allows the separation of different types of chemicals
- Locking system



TECHNICAL FEATURES

Model	ASVM11	Shipping dimensions HxWxD (mm)	2025 x 1200 x 800
Cabinet materials	10/10th epoxy painted steel structure	Ventilation	Inlet located at the bottom of the cabinet, outlet located at the top
Doors color	white RAL 9010	Openings diameter (mm)	100
Number of doors	2	Number of compartments	2
Doors type	swing doors	Containment sump	2
Hinges	right & left		115.00
Door closure	manual closing	HS Code	94032080

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Retention shelf (Réf : EE2)
- Containment sump in PVC (Réf : SBPVC)
- Shelf protection in PVC (Réf : SEPVC)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
ASVM11	1825 x 1100 x 545	2x (1800 x 540 x 475)	260		6	50	115.00