

2- glass door tall safety cabinet - 4 compartments

COMPLIANCE

ADVANTAGES

- Made in France
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Allows the separation of different types of chemicals
- Glass doors for an immediate view of the products
- Locking system
- PVC protections for acid compartment (included)



TECHNICAL FEATURES

Model	ASVM22	Ventilation	Inlet located at the bottom of the cabinet, outlet located at the top 100
Cabinet materials	10/10th epoxy painted steel structure	Openings diameter (mm)	
Doors color	white RAL 9010	Number of compartments	4
Number of doors	2	Shelves in PVC	2
Doors type	swing doors	Containment sump	3
Hinges	right & left	Containment sump in PVC	1
Door closure	manual closing	HS Code	130.00
Shipping dimensions HxWxD (mm)	2025 x 1200 x 800		94032080

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Retention shelf (Réf : EE2)
- Containment sump in PVC (Réf : SBPVC)
- Shelf protection in PVC (Réf : SEPVC)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
ASVM22	1825 x 1100 x 545	2x (1800 x 540 x 475)	220		8	50	130.00