

2- door tall safety cabinet for flammables

COMPLIANCE


FM - OSHA - NFPA

ADVANTAGES

- NFPA and OSHA compliant
- FM compliant
- Storage of flammables

BUT ALSO...

- Continious piano hinges
- Containment sump at the bottom
- NFPA, OSHA and FM compliant
- Rugged construction with double wall
- Adjustable shelves in height
- Holes for the connection to the ventilation
- Integrated door handle
- Fixation for earthing



TECHNICAL FEATURES

| | | | |
|-------------------|------------------------------------|--------------------------------|------------------------------|
| Model | MC90 | Shipping dimensions HxWxD (mm) | 1850 x 1200 x 800 |
| Compliance | FM - OSHA - NFPA | Ventilation | Both inlet and outlet |
| Cabinet materials | Double wall steel structure | Openings diameter (mm) | 50 |
| Doors color | yellow | Number of compartments | 1 |
| Number of doors | 2 | Containment sump | 1 |
| Doors type | swing doors | HS Code | 177.00 |
| Hinges | right & left | | 2147483647 |
| Door closure | manual closing | | |

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Shelf for cabinets MC90 - SC90 (Réf : ES90)

| Model | External dimensions HxWxD (mm) | Internal dimensions HxWxD (mm) | Storage capacity (L) | Containment sump | Shelves | Load capacity | Weight (kg) |
|-------|--------------------------------|--------------------------------|----------------------|------------------|---------|---------------|-------------|
| MC90 | 1651 x 1092 x 864 | 1537 x 1003 x 771 | 340 | | 2 | 159 | 177.00 |