

2- door tall safety cabinet Type 60 - 4 compartments

COMPLIANCE

EN 14470-1, EN 1363-1

ADVANTAGES

- 60-minute fire resistant and made in France
- EN 14470-1 and EN 1363-1 compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- EN 14470-1 Type 60 compliant
- Built in cellular PVC for acids-bases-toxics compartments
- Fixation for earthing (for flammable part)
- Self-closing door if the temperature rises to 47°C
- Holes for the connection to the ventilation



TECHNICAL FEATURES

Model	MR6+	Shipping dimensions HxWxD (mm)	2400 x 1300 x 800
Compliance	EN 14470-1, EN 1363-1	Ventilation	Both inlet and outlet
Cabinet materials	Corrosive part in cellular PVC, Flammable part in white RAL 9010	Openings diameter (mm)	100
Doors color	4	Number of compartments	4
Number of doors	swing doors	Containment sump	1
Doors type	right & left	HS Code	390.00
Hinges	manual & self-closing		94032080
Door closure			

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Containment sump built in polypropylene (Réf : E8A2)
- Containment sump for 1- door cabinets Types 30, 60 & 90 (Réf : B148)
- Perforated cover plate for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : C148)
- Retention shelf for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : E48)
- Sliding drawer for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : T148)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
MR6+	2205 x 1270 x 615		210		3	30	390.00