



2- door working cover safety cabinet

COMPLIANCE

NF EN ISO 1182

ADVANTAGES

- 2- point locking system
- NF EN ISO 1182 compliant
- Storage of flammables and toxics

BUT ALSO...

- Non-flammable insulation between the walls
- Normalized labelling provided with the equipment
- Holes for the connection to the ventilation
- 2- point locking system
- Containment sump at the bottom
- Self-closing doors



TECHNICAL FEATURES

Model	S1204T	Shipping dimensions HxWxD (mm)	1500 x 1200 x 800
Compliance	NF EN ISO 1182	Ventilation	Both inlet and outlet
Cabinet materials	12/10th epoxy painted steel structure	Openings diameter (mm)	100
Doors color	white RAL 9010	Number of compartments	1
Number of doors	2	Containment sump	1
Doors type	swing doors	HS Code	99.00
Hinges	right & left		94032080
Door closure	automatic closing		

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Containment sump in PVC (Réf : BP04)
- Shelf protection in PVC (Réf : EP04)
- Additional retention shelf (Réf : G04)
- Additional containment sump (Réf : GL04)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
S1204T	1295 x 1000 x 450	1145 x 950 x 400	150		3	75	99.00