

## 2- door underbench safety cabinet for flammables

### COMPLIANCE

FM - OSHA - NFPA

### ADVANTAGES

- NFPA and OSHA compliant
- FM compliant
- Storage of flammables

### BUT ALSO...

- Continuous piano hinges
- Containment sump at the bottom
- NFPA, OSHA and FM compliant
- Rugged construction with double wall
- Adjustable shelves in height
- Holes for the connection to the ventilation
- Integrated door handle
- Fixation for earthing
- Self-closing doors



### TECHNICAL FEATURES

Model	<b>SB18</b>	Shipping dimensions HxWxD (mm)	<b>810 x 1200 x 800</b>
Compliance	<b>FM - OSHA - NFPA</b>	Ventilation	<b>Both inlet and outlet</b>
Cabinet materials	<b>Double wall steel structure</b>	Openings diameter (mm)	<b>50</b>
Doors color	<b>yellow</b>	Number of compartments	<b>1</b>
Number of doors	<b>2</b>	Containment sump	<b>1</b>
Doors type	<b>swing doors</b>	HS Code	<b>58.00</b>
Hinges	<b>right &amp; left</b>		<b>2147483647</b>
Door closure	<b>self-closing</b>		

### ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

### OPTIONS

- Shelf for cabinets MB18 - SB18 (Réf : ES18)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
SB18	610 x 1092 x 457	515 x 1005 x 390	64		1	159	58.00